



[ SECURE. FAST. EFFECTIVE. ]

## V-Loc™ Absorbable Wound Closure Device Product Codes and Sizes

The V-Loc™ absorbable wound closure device is a revolutionary technology that supports optimal patient outcomes by closing wounds securely without the need to tie knots.

This guide showcases the available sizes with all applicable needles.



COVIDIEN

## V-LOC™ 90

Absorbable Wound Closure Device



PREMIUM REVERSE CUTTING ▼	Length	Color	Units Per Box	SIZE	5-0	4-0	3-0	2-0	0
				≤ 1	1	1.5	2	3	3.5
3/8 Circle 19mm  P-12	6" 15 cm	Undyed	12		VLOCM0003 ❖	VLOCM0004 ❖	VLOCM0005		
	12" 30 cm	Undyed	12		VLOCM0013 ❖	VLOCM0014 ❖	VLOCM0015		
	18" 45 cm	Undyed	12		VLOCM0023 ❖	VLOCM0024 ❖	VLOCM0025		
	24" 60 cm	Undyed	12		VLOCM0033 ❖	VLOCM0034 ❖	VLOCM0035		
3/8 Circle 24mm  P-14	12" 30 cm	Undyed	12		VLOCM0113 ❖	VLOCM0114	VLOCM0115		
	18" 45 cm	Undyed	12		VLOCM0123 ❖	VLOCM0124	VLOCM0125		
	24" 60 cm	Undyed	12		VLOCM0133 ❖	VLOCM0134	VLOCM0135		
TAPER POINT ⊙									
3/8 Circle 17mm  CV-15	6" 15 cm	Violet	12		VLOCM1203	VLOCM1204 ❖			
1/2 Circle 17mm  CV-23	12" 30 cm	Violet	12		VLOCM0813 ❖	VLOCM0814			
	18" 45 cm	Violet	12		VLOCM0823 ❖	VLOCM0824			
1/2 Circle 26mm  V-20	6" 15 cm	Violet	12		VLOCM0603 ❖	VLOCM0604			
	12" 30 cm	Violet	12		VLOCM0613 ❖	VLOCM0614			
	18" 45 cm	Violet	12		VLOCM0623 ❖	VLOCM0624	VLOCM0625		



COVIDIEN

# V-LOC™ 180

Absorbable Wound Closure Device



PREMIUM REVERSE CUTTING ▼	Length	Color	Units Per Box	SIZE	5-0	4-0	3-0	2-0	0
				≤ 1	1	1.5	2	3	3.5
3/8 Circle 19mm 	6" 15 cm	Clear	12		VLOCL0003 ♦	VLOCL0004			
	12" 30 cm	Clear	12		VLOCL0013 ♦	VLOCL0014	VLOCL0015		
	18" 45 cm	Clear	12		VLOCL0023 ♦	VLOCL0024	VLOCL0025		
	24" 60 cm	Clear	12		VLOCL0033 ♦	VLOCL0034	VLOCL0035		
3/8 Circle 24mm 	12" 30 cm	Clear	12		VLOCL0113 ♦	VLOCL0114	VLOCL0115		
	18" 45 cm	Clear	12		VLOCL0123 ♦	VLOCL0124	VLOCL0125		
	24" 60 cm	Clear	12		VLOCL0133 ♦	VLOCL0134	VLOCL0135		
TAPER POINT ◎									
1/2 Circle 17mm 	6" 15 cm	Green	12		VLOCL0803 ♦	VLOCL0804			
	12" 30 cm	Green	12		VLOCL0813 ♦	VLOCL0814			
	18" 45 cm	Green	12			VLOCL0824			
1/2 Circle 26mm 	6" 15 cm	Green	12		VLOCL0603	VLOCL0604 ♦			
	12" 30 cm	Green	12		VLOCL0613	VLOCL0614 ♦	VLOCL0615		
	18" 45 cm	Green	12			VLOCL0624 ♦	VLOCL0625		
1/2 Circle 37mm 	12" 30 cm	Green	12				VLOCL0315	VLOCL0316	
	18" 45 cm	Green	12				VLOCL0325	VLOCL0326	

	V-Loc™ 90 device	V-Loc™ 180 device
Tensile Strength	7 days 90%, 14 days 75%	7 days, 80%; 14 days, 75%; 21 days, 65%
Absorption Profile	90-110 days	180 days
Procedural Applications	Soft tissue approximation where support is required consistent with the absorption profile	Soft tissue approximation where support is required consistent with the absorption profile
Color	Undyed, violet 	Clear, green 
Composition	Glycolide, dioxanone and trimethylene carbonate	Copolymer of glycolic acid and trimethylene carbonate
Sterilization	Ethylene Oxide	Ethylene Oxide
Indications	V-Loc™ 90 device and V-Loc™ 180 absorbable wound closure devices are indicated for soft tissue approximation where use of an absorbable suture is appropriate. The safety and effectiveness of the V-Loc™ 90 device and V-Loc™ 180 device has not been established for use in fascial closures (abdominal wall, thoracic, extremity fascial closures), gastrointestinal anastomoses, cardiovascular anastomoses, neurological, ophthalmic, orthopedic or microsurgery. (Please refer to the individual package IFUs for complete instructions, indications, contraindications, warnings and precautions.)	



COVIDIEN, COVIDIEN with logo, Covidien logo and "positive results for life" are U.S. and/or internationally registered trademarks of Covidien AG.  
All other brands are trademarks of a Covidien company. © 2010 Covidien. All rights reserved. - M100533 EMEA - 04/2010

COVIDIEN (UK)  
COMMERCIAL LTD

154 FAREHAM ROAD  
GOSPORT HAMPSHIRE  
PO 13 OAS  
UNITED KINGDOM (UK)

+44 (0) 1329 224226 [T]  
+44 (0) 1329 224334 [F]

WWW.COVIDIEN.COM